

LGX-50

LIQUID GAS EXTRACTOR

The most effective way to automate your liquid gas analysis



Advantages

- No more manual prep
- Automatic internal standard addition - liquid or gas
- Sample is never exposed to the atmosphere
- System can also be used for VOC screening

Overall Value


- Best value on the market
- Proudly made in the USA

Buy With Confidence


- EST Analytical is a worldwide leader in GC sample introduction

EST  *Rely On Us*
analytical

Methane Curve



Gas Standard Overlay



Advantages

The saturated gas curves were linear with results of 0.999 or better linear regression. Carryover after the high standard was less than 0.10% and the precision of the IS and the gases tested was better than 10% RSD.



EST 
analytical *Rely On Us*

Visit
estanalytical.com